


RSAD2 Polyclonal Antibody

Catalog No.	E-AB-65580	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of U2OS cells using RSAD2 Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human RSAD2 (NP_542388.2).
GeneID	91543
Swissprot	Q8WXG1
Synonyms	RSAD2,2510004L01Rik,cig33,cig5,vig1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

RSAD2 (Radical S-Adenosyl Methionine Domain Containing 2) is a Protein Coding gene. Diseases associated with RSAD2 include Yellow Fever and Chikungunya. Among its related pathways are Interferon gamma signaling and Innate Immune System. Gene Ontology (GO) annotations related to this gene include protein self-association and iron-sulfur cluster binding.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.