


UBE2U Polyclonal Antibody


Catalog No.	E-AB-62576	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human liver using UBE2U Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse brain using UBE2U Polyclonal Antibody at dilution of 1:100 (40x lens).

Immunogen Information

Immunogen	Recombinant fusion protein of human UBE2U (NP_689702.1).
GeneID	148581
Swissprot	Q5VVX9
Synonyms	UBE2U

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:100

Background

UBE2U (Ubiquitin Conjugating Enzyme E2 U) is a Protein Coding gene. Diseases associated with UBE2U include Brachydactyly, Type E2. Among its related pathways are Remodeling of Adherens Junctions and Innate Immune System. Gene Ontology (GO) annotations related to this gene include ligase activity and ubiquitin protein ligase activity. An important paralog of this gene is UBE2B.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.