


PSIP1 Polyclonal Antibody

Catalog No.	E-AB-62378	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using PSIP1 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human PSIP1 (NP_150091.2).
GeneID	11168
Swissprot	O75475
Synonyms	PSIP1,DFS70,LEDGF,PAIP,PSIP2,p52,p75

Product Information

Calculated MW	37kDa/60kDa
Observed MW	90kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

PSIP1 (PC4 And SFRS1 Interacting Protein 1) is a Protein Coding gene. Diseases associated with PSIP1 include Dermatitis, Atopic and Immune Deficiency Disease. Among its related pathways are Neuroscience and NF-KappaB Family Pathway. An important paralog of this gene is HDGFL2.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.