


PTCD3 Polyclonal Antibody

Catalog No.	E-AB-62466	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using PTCD3 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human PTCD3 (NP_060422.4).
GeneID	55037
Swissprot	Q96EY7
Synonyms	PTCD3,MRP-S39

Product Information

Calculated MW	31kDa/78kDa
Observed MW	79kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

PTCD3 (Pentatricopeptide Repeat Domain 3) is a Protein Coding gene. Diseases associated with PTCD3 include Leigh Syndrome, French Canadian Type and Leigh Syndrome. Among its related pathways are Organelle biogenesis and maintenance and Mitochondrial translation. Gene Ontology (GO) annotations related to this gene include ribosomal small subunit binding.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.