


ARL8B Polyclonal Antibody

Catalog No.	E-AB-65602	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Confocal immunofluorescence analysis of U2OS cells using ARL8B Polyclonal Antibody at dilution of 1:100. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human ARL8B (NP_060654.1).
GeneID	55207
Swissprot	Q9NVJ2
Synonyms	ARL8B,ARL10C,Gie1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

ARL8B (ADP Ribosylation Factor Like GTPase 8B) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include GTP binding and obsolete signal transducer activity. An important paralog of this gene is ARL8A.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.