


EBNA1BP2 Polyclonal Antibody

Catalog No.	E-AB-63939	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Confocal immunofluorescence analysis of U2OS cells using EBNA1BP2 Polyclonal Antibody at dilution of 1:200. Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human EBNA1BP2 (NP_006815.2).
GeneID	10969
Swissprot	Q99848
Synonyms	EBNA1BP2,EBP2,NOBP,P40

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:200

Background

EBNA1BP2 (EBNA1 Binding Protein 2) is a Protein Coding gene. Diseases associated with EBNA1BP2 include Robinow Syndrome, Autosomal Dominant 1 and Intracranial Hypotension.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.