


ZRANB3 Polyclonal Antibody

Catalog No.	E-AB-65632	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using ZRANB3 Polyclonal Antibody at dilution of 1:3000.

Immunogen Information

Immunogen	Recombinant fusion protein of human ZRANB3 (NP_001273497.1).
GeneID	84083
Swissprot	Q5FWF4
Synonyms	ZRANB3,4933425L19Rik,AH2

Product Information

Calculated MW	42kDa/49kDa/67kDa/123kDa
Observed MW	123kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:1000-1:3000

Background

ZRANB3 (Zinc Finger RANBP2-Type Containing 3) is a Protein Coding gene. Diseases associated with ZRANB3 include Warburg Micro Syndrome and Warburg Micro Syndrome 1. Gene Ontology (GO) annotations related to this gene include nucleic acid binding and endonuclease activity. An important paralog of this gene is CHD7.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.