


ROBO4 Polyclonal Antibody

Catalog No.	E-AB-62460	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using ROBO4 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human ROBO4 (NP_061928.4).
GeneID	54538
Swissprot	Q8WZ75
Synonyms	ROBO4,ECSM4,MRB

Product Information

Calculated MW	75kDa/98kDa/107kDa
Observed MW	130kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

ROBO4 (Roundabout Guidance Receptor 4) is a Protein Coding gene. Diseases associated with ROBO4 include Aortic Valve Disease 3 and Diaphragmatic Hernia, Congenital. Among its related pathways are Development Slit-Robo signaling. An important paralog of this gene is SDK1.

For Research Use Only

Thank you for your recent purchase.
If you would like to learn more about antibodies, please visit www.elabscience.com.

Focus on your research Service for life science

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.