


UBE2Q2 Polyclonal Antibody

Catalog No.	E-AB-62558	Reactivity	H,M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using UBE2Q2 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human UBE2Q2 (NP_775740.1).
GeneID	92912
Swissprot	Q8WVN8
Synonyms	UBE2Q2

Product Information

Calculated MW	36kDa/38kDa/40kDa/42kDa
Observed MW	43kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

UBE2Q2 (Ubiquitin Conjugating Enzyme E2 Q2) is a Protein Coding gene. Diseases associated with UBE2Q2 include Hypophosphatemic Nephrolithiasis/Osteoporosis. Among its related pathways are Metabolism of proteins and Innate Immune System. Gene Ontology (GO) annotations related to this gene include ligase activity and acid-amino acid ligase activity. An important paralog of this gene is UBE2Q1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications: WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.