


NEDD8 Polyclonal Antibody


Catalog No.	E-AB-40240	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human liver using NEDD8 Polyclonal Antibody at dilution of 1:200



Immunohistochemistry of paraffin-embedded Human cerebral cortex using NEDD8 Polyclonal Antibody at dilution of 1:200(Elabscience® Product Detected by Lifespan).

Immunogen Information

Immunogen	Recombinant Human NEDD8 protien
Swissprot	Q15843
Synonyms	FLJ43224,MGC104393,MGC125896,MGC125897,NED8,Nedd8,Neddylin,Rub1,Ubiquitin like protein Nedd 8

Product Information

Buffer	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.
Purify	Affinity purification
Dilution	IHC 1:50-1:200

Background

NEDD8 (neural precursor cell expressed, developmentally down-regulated 8) is a protein-coding gene. Diseases associated with NEDD8 include glomuvenous malformation. GO annotations related to this gene include ubiquitin protein ligase binding.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.