


beta Tubulin Polyclonal Antibody

Catalog No.	E-AB-40483	Reactivity	M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western Blot analysis of Mouse intestine, Mouse kidney and Rat testis using beta Tubulin Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Synthetic peptide of Mouse Tubulin beta-5 chain
GeneID	22154
Swissprot	P99024
Synonyms	TUBB, TUBB1, TUBB2, TUBB5

Product Information

Calculated MW	50 kDa
Observed MW	55 kDa
Buffer	PBS with 0.05% Proclin300, 50% glycerol, pH7.3.
Purify	Antigen Affinity Purification
Dilution	WB 1:500-1:2000

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.